

The Product

CHEM R-OL is a high grade litho IPA Substitute which is added in the fountain solution mixture.

Features

- Helps in optimisation of print density and gloss.
- Ensures optimum dosage and sharper dot reproduction.
- Eliminates the use of IPA in printing.
- Low VOC content.
- Compatible with all types of alcohol dampening systems.

User Guidelines

- Always ensure that CHEM R-OL is added to the mixed fount solution and not to the concentrated fount.
- Dosage to be at 2% - 3% depending on the press condition.
- Continuous replenishment is necessary to reduce evaporation loss and maintain dosage.
- Recommended to be temperature controlled.

Safety & Security

- Always make use of Protective gloves and goggles (eye wear) while using the wash.
- Keep away from naked flame or dangerous flammable and hazardous material.
- Store at moderate temperatures in a dry, well-ventilated area.
- Avoid direct skin and eye contact. Do not inhale the vapours / fumes directly.
- In case of Skin / Eye contact immediately rinse with water thoroughly and consult a physician.
- In case of ingestion immediately consult a physician. - Do not induce vomiting.

Available In

- | | | |
|--------------|--------------|-----------------------------|
| ● 20 / 25 L | ● 5 L | ● 1000 L |
| Plastic Cont | Plastic Cont | (Packing only avail on req) |

FOUNTAIN SOLUTION ADDITIVES